




Prairie Ceramics

HIGH TEMPERATURE FURNACE SOLUTIONS

Standard High Purity Oxide TGA Crucibles



Prairie Ceramics provides high performance standard form crucibles for use in Thermal-gravimetric Analysis (TGA) machines. Sizes conform to commonly used shapes and are available in high purity 99.9% alumina and Prairie's exclusive Yttria-Super (99.9% Pure) body.

See standard sizes available below or contact us for a specific size required in your process.

P/N (in mm)	OD	HT	WALL	LT	HD
TGA◇-050-080-**	5.0	8.0	.7	1	1
TGA◇-055-040-**	5.5	4.0	.7	1	1
TGA◇-060-020-**	6.0	2.0	.7	1	1
TGA◇-068-040-**	6.8	4.0	.7	1	1
TGA◇-106-030-**	10.6	3.0	.7	1	1

◇ Use "C" for Crucible Only, "S" for Crucible and Lid Set

** Use YS for Yttria, AD999 for 99.9% alumina, YTZP for Zirconia

